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

Canonical mechanisms of GPCR signalling

a G protein-mediated signalling by GPCRs

b Arrestin-mediated signalling by GPCRs

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


- Bioluminescent sensor was genetically modified form of luciferase .
- Its luminescence activity dependent on reversible allosteric interaction with ligand.
- It can monitor cAMP real time in live-cells kinetically.

WorkFlow

Agonist

Antagonist


Read Luminescence to Determine Post Incubation Time Signal


Read Luminescence to Determine Post Incubation Time Signal

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

β -Adrenergic Receptor Agonist Titration

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


分析	靶点
AlamarBlue®	Cell Proliferation, Viability & Counting
Amplex® Red	Enzyme Activity
phosphoELISA	Kinases
Vybrant assay C12 Rezarurin	Metabolic Labeling
Vybrant MTT Cell Proliferation Assay Kit	Cell Counting
Lipid Peroxidation Assay Kit (using MDA)	Oxidative Stress
CellTiter 96® Cell Proliferation Assay (MTT)	Cell Proliferation
DC Protein, Bradford, Lowry & BCA Protein Assays	Protein Quantification
OxiSelect™ Superoxide Dismutase Activity Assay	Oxidative Stress
Glycoprotein Carbohydrate Estimation Kit	Glycoprotein Quantification
IDEXX FlockChek™ Avian Influenza Multi-species ELISA Assay	Avian Influenza
Melamine ELISA Assay	Food Safety
CaspACE™ Assay	Apoptosis
Growth Assay (OD at 600nm)	Cell Growth
Haemoglobin denaturation Assay	Haemoglobin

分析	靶点
Adapta [®] Universal Kinase Assay	Kinases
AlamarBlue [®]	Cell Proliferation, Viability & Counting
Amplex [®] Red	Enzyme Activity
EnzCheck	Proteases
Fluo-4 AM	GPCR
LIVE/DEAD [®] Viability/Cytotoxicity Kit	Apoptosis
PiPer [™] Phosphate Assay Kit	Phosphatases
SYTOX [®] stains	Viability
Vivid [®] CYP450	Cytochrome P450
Quant-iT [™] PicoGreen [®] dsDNA Assay Kit	Nucleic Acid Quantification
Transcreener [®] GDP FI Assay	GTPase
FluxOR [™] Potassium Ion Channel Assay Kits	Ion Channels
OxiSelect [™] ROS Assay Kit	Oxidative Stress
CellTiter-Fluor [™] Cell Viability Assay	Cell Viability
FluoroBlok Cell Invasion Assays (using Calcein AM)	Cell Migration

分析	靶点
CellTiter-Glo [®] Luminescent Cell	Cell Proliferation
Dual-Luciferase [®] Reporter (DLR [™]) Assay	Gene Expression
P450-Glo [™] Assay	CYP Activities
ADP-Glo [™] Assay	Kinase
Caspase-Glo [®] 3/7 Assay	Apoptosis
cAMP-Glo [™] Assay	GPCR
GloSensor [™] cAMP Assay	GPCR
Aequorin	GPCR
ATPlite [™] Luminescence Assay	Cell Proliferation/Kinases
aCella -TOX [™] Bioluminescence Cytotoxicity Assay	Cell Death
ToxiLight [®] Non-destructive Cytotoxicity BioAssay Kit	Cell Death/Cell Proliferation
ATP quantification assay (luciferine-luciferase)	ATP Quantification
Britelite [™] plus Reporter Gene Assay	Gene Expression
Steadylite plus [™] Reporter Gene Assay	Gene Expression
Neolite Reporter Gene Assay	Gene Expression

	TRF	
	LANCE/HTRF	DELFI A
GPCR 功能分析:		
cAMP		✓
IP ₃		
GTP Binding		✓
受体结合分析:		
GPCRs	✓	✓
Nuclear	✓	✓
酶分析:		
Tyrosine Kinase	✓	✓
Serine/Threonine Kinase	✓	✓
Phosphatase	✓	
Protease	✓	
Transferase	✓	✓
Polymerase	✓	✓
Phospholipase		
Helicase	✓	
Phosphodiesterase		
蛋白相互作用:		
Protein:Protein	✓	✓
Protein:Peptide	✓	✓
Protein:DNA		✓
其他功能分析:		
Cell Proliferation/Viability		✓

分析	靶点
PolarScreen™	Kinases
PolarScreen™ Competitive Binding Assays	Nuclear Receptors
Transcreener® AMP/GMP Assay	PDE's
Transcreener® GDP FP Assay	GTPase
Transcreener® UDP Assay	UDP-glycosyltransferases
Transcreener® ADP2 FP Assay	Kinases
[FP] ²	cAMP

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Memorial
University of Newfoundland

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In-vitro cell-based assays that utilize a 3D structure of aggregated cells using hydrogels, scaffolds or spheroids to provide more physiologically relevant data.

ADME/Tox (21)

In-vitro assays using e.g. microplates for drug absorption, drug metabolism

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Microplate assays used for biologic drug discovery and development: clonal selection, characterization, formulation studies and immunogenicity.

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Demonstration of automated assays

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A wide range of cell-based assays not in the context of drug discovery.

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Assays using purified components

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Fluorescence microscopy of cells and phenotypic screening using microplates.

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Assays assisting in R&D of biofuels

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Cell-based, functional assays used typically in Hit ID and Hit to Lead small molecule screening efforts.

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Quantification of biological markers in biofluids (i.e. serum/plasma)

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Kinetic studies using purified enzymes in a microplate format.

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Cell-based and biochemical assays to measure activity and small molecule modulation of epigenetic target proteins involved with post-translational modifications of histone tails.

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Assays measuring contaminants in food and beverages including bacteria and processing by-products that influence customer health and product quality.

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